

# **Premier Wash Reagent**



01-03-0098

Fill 940mL

Pour d'autres langues Für andere Sprachen Para otras lenguas Per le altre lingue Dla innych języków

Para outras línguas Για τις άλλες λώσσες

För andra språk For andre språk



#### INTENDED USE

 $\underline{\text{Hemoglobin A1c}} - \text{This reagent is intended for use with the Trinity Biotech Premier Hb9210}^{\text{TM}} \ \text{HbA1c}$ 

The Premier system is intended for the quantitative measurement of hemoglobin A1c (HbA1c) in human capillary and venous whole blood. HbA1c is used for the monitoring of long-term glycemic control in individuals with diabetes mellitus. For in vitro diagnostic use only.

## SUMMARY AND EXPLANATION OF TEST

Hemoglobin A1c - Assessment of hemoglobin A1c has proven useful in the control of diabetes. Reagents are performance validated to assure accuracy and precision with the Trinity Biotech assay and system for the measurement of hemoglobin A1c.

Reagent is ready for use. No substitutions or other uses are permitted, registered or authorized.

## STORAGE AND STABILITY



Store at ambient temperatures (2 – 28°C). Do not allow to freeze.

The reagent is stable until the expiration date indicated on the label when kept tightly closed and protected from extreme environmental conditions.




See the container label for the expiration date. DO NOT USE after the expiration date.

#### PRECAUTIONS

For in vitro diagnostic use only. Avoid skin contact. Consult the product SDS for safety information. This reagent is used in conjunction with blood testing equipment and warrants handling under universal precaution procedures for safety.

ORDERING INFORMATION		
Reference No.	Item	Quantity
01-03-0098	Premier Wash	940mL





Trinity Biotech Kansas City, MO 64132 USA Tel. 1 800-325-3424 Fax: 1816-361-1974



Trinity Biotech plc Bray Co. Wicklow, Ireland Tel. 353 1 2769800 Fax 353 1 2769888 www.trinitybiotech.com

